




SURGIMEDIC

Quality. Commitment. Innovation.



Complete Wound Closure Catalog

About us ...

When you depend upon the quality and effectiveness of the suture in your hand, **SURGIMEDIC** is the names you can trust. Our proprietary manufacturing process gives us complete control over the quality of every product that we sell. Today, **SURGIMEDIC** brand sutures continue to meet the demand of the newest techniques in office based and dental surgery procedures.

For over ten years, the **SURGIMEDIC** name has been synonymous with providing the highest quality needles and sutures to the dermatologic, plastic, microsurgical, and ophthalmic communities. **SURGIMEDIC** isn't just a brand name; it's our reputation. Our goal is to provide innovative, high quality products and service that will ensure the continued success of the surgical professionals we serve.



Quality.

Superior design, materials and manufacturing deliver the confidence you expect from the products you need.

SURGIMEDIC EXPERT tools and expert team provide all the required support when you need it, making us easy to partner with.



Commitment.

We listen, assess and constantly seek to improve our products and services to support you in delivering your patient and management objectives.

Our **SURGIMEDIC** team is expert and we have a **SURGIMEDIC EXPERT** tool kit to match, delivering best value today and tomorrow.




Innovation.

Unique and speciality products deliver the performance and convenience you need.

SURGIMEDIC EXPERT provides the accurate information you need, assisting you with clinical choice and management efficiency.





SUTURE TYPE **CODE**

SUTCRYL POLYGLACTIN 910	T
SURCRYL POLYCLYCOLIC ACID	G
SURCARPONE POLYGLICARPONE 25	D
SURCANONE POLYDIOXANONE	C
SURLENE POLYPROPYLENE (monofilament)	P
SURLON NYLON (monofilament)	N
SURSILK SILK BRAIDED	S
SURLYEST POLYESTER BRAIDED	L

CROSS SECTION **CODE**

● round bodied	R
▼ reverse cutting	C
▲ cutting	G
⊗ tabercut	T
◄ micro-point spatula curved	P
▲ straight cutting	S

SUTURE LENGTH **CODE**

15 cm	015
30 cm	030
45 cm	045
60 cm	060
75 cm	075
90 cm	090
100 cm	100
150 cm	150
180 cm	180

SUTURE DIAMETER **CODE**



8/0	80
7/0	70
6/0	60
5/0	50
4/0	40
3/0	30
2/0	20
0	00
1	01
2	02

NEEDLE LENGTH **CODE**






13	13
15	15
16	16
18	18
19	19
21	21
25	25
27	27
30	30
33	33
35	35
37	37
40	40
43	43
45	45
48	48
50	50

NEEDLES



SURGIMEDIC's eyeless needles are made from special stainless steel with optimal flexibility and strength. Special surface treatment and precision grinding of the point or edge ensure minimal resistance on insertion and easy passage of the needle through the tissue. The firmness with which the needle is attached to the suture is tested in accordance with the regulations of harmonized standards for surgical suturing materials according to the European Pharmacopoeia. We recommend use of needles with welded suture to minimize wound trauma due to double thread.





- **ROUND BODIED**
Designed to provide easy penetration of intended tissue. Widely used in soft tissue where separation of tissue fibers is preferred, rather than cutting the fibers apart.
- ▲ **CURVED CUTTING**
Characterized by its triangular cutting pattern, with the apex cutting edge on the inside of the needle. Widely used in reconstructive plastic surgery and in cuticular closing.
- ▼ **REVERSE CUTTING**
Characterized by its triangular cutting pattern, with the apex cutting edge on the outside of the needle curvature. Widely used in reconstructive plastic surgery and in cuticular closing.
- ⊕ **TAPERCUT**
A round bodied needle with a cutting point; this needle gives good and easy penetration in tough tissues with little tissue trauma. Widely used for tendons, orthopedics and sternum closure.
- ▼ **MICRO-POINT SPATULA CURVED**
This type of needle is a special needle for ophthalmic surgery. the needle is flat, thin and has an extremely sharp cutting edge, which makes it suitable for penetration between the layers of scleral or corneal tissue.
- ⊙ **BLUNT TAPER POINT**
Cylindrical bodied needle with a rounded point. Particularly designed for suturing soft tissue such as in liver and kidneys.




Absorbable Sutures



Non-Absorbable Sutures




Absorption Times




1. Brand name
2. Product material
3. EP size
4. USP size
5. Thread length
6. Product code
7. Sterile product using Ethylene Oxide
8. Manufacture Name and address
9. Authorised representative in the European Community
10. Do not re-sterilize
11. Caution symbol
12. Keep away from sunlight symbol
13. Keep dry symbol
14. Upper limit of storage temperature (25 °C)
15. Do not re-use
16. CE mark and number
17. Country of origin
18. Expiration date
19. Manufacturing date
20. LOT number
21. Needle drawing
22. Needle curvature
23. Needle length
24. Needle point drawing
25. Needle point type
26. Number of sutures


1. Suture (Inner Pouch) Label



2. Inner Carton (1 Dozen) Label



3. Shipper Carton (96 Dozens) Label



Absorbable Polyglycolic Acid Sutures



Type
Braided and coated multifilament

Composition
Polyglycolic acid

Coating
Polycaprolactone and calcium stearate (1%)

Color
Beige - undyed

Tissue reaction
Minimal

Absorption
The hydrolytic action by which the material is broken down results in total absorption in approximately 42 days. Approximately 50% of tensile strength remain after one week.

Distinctive characteristics
Fast absorbing polyglycolic acid is intended to model the performance of collagen (surgical gut) suture.
High tensile strength.
Excellent handling.
Excellent knot security.
Excellent tissue support.
Smooth passage through the tissues.


Indications
Superficial soft tissue approximation of the skin and mucosa, where only short-term wound support is required.

especially in
Plastic surgery
Obstetrics (Episiotomy repair)
Gynaecology
Urology (Circumcision)
Ophthalmic surgery (Conjunctiva)
Pediatric surgery

Contraindication
Due to rapid loss of tensile strength, this suture should not be used where extended approximation of tissue under stress is required or where wound support beyond 7 days is required.








Packaging
SURCRYL rapid is available in USP sizes 6/0 through 1 (metric sizes 0.7-4.0).
SURCRYL rapid is available with needles in one dozen boxes.
Three dozen suture boxes available on customer request.






TENSILE STRENGTH



SURCRYL rapid


Absorbable Polyglycolic Acid Sutures

Round bodied 1/2 circle	Metric	USP	Length	Color	Order number
 16MM	0.7	6/0	75 cm	undyed	RR6016075
	1	5/0	75 cm	undyed	RR5016075
 19MM	1	5/0	75 cm	undyed	RR5019075
	1.5	4/0	75 cm	undyed	RR4019075
 25MM	1.5	4/0	75 cm	undyed	RR4025075
	2	3/0	75 cm	undyed	RR3025075
 30MM	2	3/0	75 cm	undyed	RR3030075
	3	2/0	75 cm	undyed	RR2030075
	3.5	0	75 cm	undyed	RR0030075
 35MM	3	2/0	75 cm	undyed	RR2035075
	3.5	0	75 cm	undyed	RR0035075
	4	1	75 cm	undyed	RR0135075
 40MM	3.5	0	75 cm	undyed	RR0040075
	4	1	75 cm	undyed	RR0140075
	5	2	75 cm	undyed	RR0240075
 50MM	3.5	0	75 cm	undyed	RR0050075
	4	1	75 cm	undyed	RR0150075
	5	2	75 cm	undyed	RR0250075

Reverse cutting 3/8 circle	Metric	USP	Length	Color	Order number
 16MM	0.7	6/0	75 cm	undyed	RC6016075
	1	5/0	75 cm	undyed	RC5016075
 19MM	1.5	4/0	75 cm	undyed	RC4019075
	2	3/0	75 cm	undyed	RC3019075
 25MM	2	3/0	75 cm	undyed	RC3025075
	3	2/0	75 cm	undyed	RC2025075
 30MM	2	3/0	75 cm	undyed	RC3030075
	3	2/0	75 cm	undyed	RC2030075
	3.5	0	75 cm	undyed	RC0030075
 35MM	3	2/0	75 cm	undyed	RC2035075
	3.5	0	75 cm	undyed	RC0035075

Absorbable Polyglycolic Acid Sutures

Cutting 1/2 circle



Metric	USP	Length	Color	Order number
1.5	4/0	45 cm	undyed	RG4019045
2	3/0	45 cm	undyed	RG3019045
3	2/0	75 cm	undyed	RG2035075
3.5	0	75 cm	undyed	RG0035075
4	1	75 cm	undyed	RG0135075
3.5	0	90 cm	undyed	RG0040090
4	1	90 cm	undyed	RG0140090
5	2	90 cm	undyed	RG0240090

Cutting straight




60MM



Metric	USP	Length	Color	Order number
2	3/0	75 cm	undyed	RS3060075

Tapercut 1/2 circle



Metric	USP	Length	Color	Order number
1.5	4/0	75 cm	undyed	RT4019075
2	3/0	75 cm	undyed	RT3019075
3	2/0	75 cm	undyed	RT2035075
3.5	0	75 cm	undyed	RT0035075
4	1	75 cm	undyed	RT0135075
3.5	0	75 cm	undyed	RT0040075
4	1	75 cm	undyed	RT0140075
5	2	75 cm	undyed	RT0240075

Without needle



Metric	USP	Length	Color	Order number
2	3/0	150 cm	undyed	RNN30150
3	2/0	150 cm	undyed	RNN20150
3.5	0	150 cm	undyed	RNN00150
4	1	150 cm	undyed	RNN01150
5	2	150 cm	undyed	RNN02150

SUTCRYL

Absorbable Polyglactine 910 Sutures



Type

Braided and coated multifilament

Composition

Polyglactine 910, a copolymer made of 90% glycolide and 10%L-lactide

Coating

Polycaprolactone and calcium stearate (1%)

Color

Violet

Tissue reaction

Minimal

Absorption

The hydrolytic action by which the material is broken down results in total absorption in approximately 56 to 70 days. The material retains approximately 75% of its initial tensile strength at the end of second week. Approximately 50% of tensile strength remain after three weeks.

Distinctive characteristics

High tensile strength.
Excellent handling.
Excellent knot security.
Excellent tissue support.
Smooth passage through the tissues.

Indications

General soft tissue closing and/or ligation.

especially in

General surgery
Skin closure
Gastrointestinal surgery
Plastic surgery
Gynaecology
Urology
Ophthalmic surgery


Contraindication

Cardiovascular surgery and neurosurgery.
Since it is an absorbable suture; It should not to use when an extended wound support is required.









Packaging



SUTCRYL is available in USP sizes 6/0 through 2 (metric sizes 0.7-5.0).
SUTCRYL is available with needles in one dozen boxes.
Three dozen suture boxes available on customer request.



TENSILE STRENGTH



Absorbable Polyglactine 910 Sutures







Round bodied 1/2 circle		Metric	USP	Length	Color	Order number
		0.7	6/0	75 cm	violet	TR6016075
		1	5/0	75 cm	violet	TR5016075
		1	5/0	75 cm	violet	TR5019075
		1.5	4/0	75 cm	violet	TR4019075
		1.5	4/0	75 cm	violet	TR4025075
		2	3/0	75 cm	violet	TR3025075
		2	3/0	75 cm	violet	TR3030075
		3	2/0	75 cm	violet	TR2030075
		3.5	0	75 cm	violet	TR0030075
		3	2/0	75 cm	violet	TR2035075
		3.5	0	90 cm	violet	TR0035090
		4	1	90 cm	violet	TR0135090
		3	2/0	75 cm	violet	TR2040075
		3.5	0	75 cm	violet	TR0040075
		4	1	90 cm	violet	TR0140090
		4	1	75 cm	violet	TR0140075
		5	2	90 cm	violet	TR0240090
		5	2	75 cm	violet	TR0240075
		3.5	0	90 cm	violet	TR0050090
		4	1	90 cm	violet	TR0150090
		5	2	90 cm	violet	TR0250090





Round bodied 5/8 circle		Metric	USP	Length	Color	Order number
		2	3/0	75 cm	violet	TR3027075

Round bodied 3/8 circle		Metric	USP	Length	Color	Order number
		2	3/0	75 cm	violet	TR3026075

SUTCRYL


Absorbable Polyglactine 910 Sutures

Reverse cutting 3/8 circle ▼	Metric	USP	Length	Color	Order number
	0.7	6/0	75 cm	violet	TC6016075
	1	5/0	75 cm	violet	TC5016075
	1.5	4/0	75 cm	violet	TC4019075
	2	3/0	75 cm	violet	TC3019075
	2	3/0	75 cm	violet	TC3025075
	3	2/0	75 cm	violet	TC2025075
	2	3/0	75 cm	violet	TC3030075
	3	2/0	75 cm	violet	TC3030075
	3.5	0	75 cm	violet	TC0030075
	3	2/0	75 cm	violet	TC2035075
	3.5	0	75 cm	violet	TC0035075
	3	2/0	75 cm	violet	TC2040075
	3.5	0	75 cm	violet	TC0040075
	4	1	75 cm	violet	TC0140075

Reverse cutting 1/2 circle ▼	Metric	USP	Length	Color	Order number
	1.5	4/0	75 cm	violet	TC4018075
	2	3/0	75 cm	violet	TC3018075
	3	2/0	75 cm	violet	TC2037075
	3.5	0	90 cm	violet	TC0037090
	4	1	90 cm	violet	TC0137090
	3	2/0	90 cm	violet	TC2045090
	3.5	0	90 cm	violet	TC0045090
	4	1	90 cm	violet	TC0145090
	5	2	90 cm	violet	TC0245090
	3.5	0	90 cm	violet	TC0048090
	4	1	90 cm	violet	TC0148090
	5	2	90 cm	violet	TC0248090

Absorbable Polyglactine 910 Sutures

Tapercut 1/2 circle



Metric	USP	Length	Color	Order number
1.5	4/0	75 cm	violet	TT4025075
2	3/0	75 cm	violet	TT3025075
3	2/0	75 cm	violet	TT2030075
3.5	0	75 cm	violet	TT0030075
3	2/0	90 cm	violet	TT2037090
3.5	0	90 cm	violet	TT0037090
4	1	90 cm	violet	TT0137090
3.5	0	90 cm	violet	TT0040090
4	1	90 cm	violet	TT0140090
5	2	90 cm	violet	TT0240090

Cutting straight



60MM



Metric	USP	Length	Color	Order number
2	3/0	75 cm	violet	TS3060075
2	2/0	75 cm	violet	TS2060075

Without needle



Metric	USP	Length	Color	Order number
2	3/0	150 cm	violet	TNN30150
3	2/0	150 cm	violet	TNN20150
3.5	0	150 cm	violet	TNN00150
4	1	150 cm	violet	TNN01150
5	2	150 cm	violet	TNN02150
2	3/0	6 X 45 cm	violet	TNN30645
3	2/0	6 X 45 cm	violet	TNN20645
3.5	0	6 X 45 cm	violet	TNN00645
4	1	6 X 45 cm	violet	TNN01645
5	2	6 X 45 cm	violet	TNN02645

SURCRYL

Absorbable Polyglycolic Acid Sutures



Type

Braided and coated multifilament

Composition

Polyglycolic acid

Coating

Polycaprolactone and calcium stearate (1%)

Color

Violet

Tissue reaction

Minimal

Absorption

The hydrolytic action by which the material is broken down results in total absorption in approximately 60 to 90 days. The material retains approximately 70% of its initial tensile strength at the end of second week. Approximately 50% of tensile strength remain after three weeks.

Distinctive characteristics

High tensile strength.
Excellent handling.
Excellent knot security.
Excellent tissue support.
Smooth passage through the tissues.

Indications

General soft tissue closing and/or ligation.

especially in

General surgery
Skin closure
Gastrointestinal surgery
Plastic surgery
Gynaecology
Urology
Ophthalmic surgery
Orthopaedics


Contraindication

Cardiovascular surgery and neurosurgery.
Since it is an absorbable suture; It should not be used when an extended wound support is required.













Packaging

SURCRYL is available in USP sizes 6/0 through 2 (metric sizes 0.7-5.0).
SURCRYL is available with needles in one dozen boxes.
Three dozen suture boxes available on customer request.

TENSILE STRENGTH













Absorbable Polyglycolic Acid Sutures

Round bodied 1/2 circle		Metric	USP	Length	Color	Order number
	16MM	0.7	6/0	75 cm	violet	GR6016075
		1	5/0	75 cm	violet	GR5016075
	19MM	1	5/0	75 cm	violet	GR5019075
		1.5	4/0	75 cm	violet	GR4019075
	25MM	1.5	4/0	75 cm	violet	GR4025075
		2	3/0	75 cm	violet	GR3025075
	30MM	2	3/0	75 cm	violet	GR3030075
		3	2/0	75 cm	violet	GR2030075
		3.5	0	75 cm	violet	GR0030075
	35MM	3	2/0	75 cm	violet	GR2035075
		3.5	0	90 cm	violet	GR0035090
		4	1	90 cm	violet	GR0135090
	40MM	3	2/0	75 cm	violet	GR2040075
		3.5	0	75 cm	violet	GR0040075
		4	1	90 cm	violet	GR0140090
		4	1	75 cm	violet	GR0140075
		5	2	90 cm	violet	GR0240090
		5	2	75 cm	violet	GR0240075
	50MM	3.5	0	90 cm	violet	GR0050090
		4	1	90 cm	violet	GR0150090
		5	2	90 cm	violet	GR0250090
Round bodied 5/8 circle		Metric	USP	Length	Color	Order number
	27MM	2	3/0	75 cm	violet	GR3027075
Round bodied 3/8 circle		Metric	USP	Length	Color	Order number
	26MM	2	3/0	75 cm	violet	GR3026075

SURCRYL


Absorbable Polyglycolic Acid Sutures

Reverse cutting 3/8 circle ▼	Metric	USP	Length	Color	Order number
	0.7	6/0	75 cm	violet	GC6016075
	1	5/0	75 cm	violet	GC5016075
	1.5	4/0	75 cm	violet	GC4019075
	2	3/0	75 cm	violet	GC3019075
	2	3/0	75 cm	violet	GC3025075
	3	2/0	75 cm	violet	GC2025075
	2	3/0	75 cm	violet	GC3030075
	3	2/0	75 cm	violet	GC2030075
	3.5	0	75 cm	violet	GC0030075
	3	2/0	75 cm	violet	GC2035075
	3.5	0	75 cm	violet	GC0035075
	3	2/0	75 cm	violet	GC2040075
	3.5	0	75 cm	violet	GC0140075
	4	1	75 cm	violet	GC0240075

Reverse cutting 1/2 circle ▼	Metric	USP	Length	Color	Order number
	1.5	4/0	75 cm	violet	GC4018075
	2	3/0	75 cm	violet	GC3018075
	3	2/0	75 cm	violet	GC2037075
	3.5	0	90 cm	violet	GC0037090
	4	1	90 cm	violet	GC0137090
	3	2/0	90 cm	violet	GC2040090
	3.5	0	90 cm	violet	GC0040090
	4	1	90 cm	violet	GC0140090
	5	2	90 cm	violet	GC0240090
	3.5	0	90 cm	violet	GC0048090
	4	1	90 cm	violet	GC0148090
	5	2	90 cm	violet	GC0248090

Absorbable Polyglycolic Acid Sutures

Tapercut 1/2 circle



Metric	USP	Length	Color	Order number
1.5	4/0	75 cm	violet	GT4025075
2	3/0	75 cm	violet	GT3025075
3	2/0	75 cm	violet	GT2030075
3.5	0	75 cm	violet	GT0030075
3	2/0	90 cm	violet	GT2035090
3.5	0	90 cm	violet	GT0035090
4	1	90 cm	violet	GT0135090
3.5	0	90 cm	violet	GT0040090
4	1	90 cm	violet	GT0140090
5	2	90 cm	violet	GT0240090

Cutting straight



60MM



Metric	USP	Length	Color	Order number
2	3/0	75 cm	violet	GS3060075
2	2/0	75 cm	violet	GS2060075

Without needle



Metric	USP	Length	Color	Order number
2	3/0	150 cm	violet	GNN30150
3	2/0	150 cm	violet	GNN20150
3.5	0	150 cm	violet	GNN00150
4	1	150 cm	violet	GNN01150
5	2	150 cm	violet	GNN02150
2	3/0	6 X 45 cm	violet	GNN30645
3	2/0	6 X 45 cm	violet	GNN20645
3.5	0	6 X 45 cm	violet	GNN00645
4	1	6 X 45 cm	violet	GNN01645
5	2	6 X 45 cm	violet	GNN02645
2	3/0	10 X 60 cm	violet	GNN301060
3	2/0	10 X 60 cm	violet	GNN201060
3.5	0	10 X 60 cm	violet	GNN001060
4	1	10 X 60 cm	violet	GNN011060

SURCAPRONE

Absorbable Polyglecaprone 25 Sutures



Type

Monofilament

Composition

Polyglecaprone 25, a copolymer made of 75% glycolide and 25% ϵ -caprolactone

Coating

None

Color

Violet

Tissue reaction

Minimal

Absorption

The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 120 days. The material retains approximately 30% of its initial tensile strength at the end of second week. Approximately 60% of tensile strength remain after one week.

Distinctive characteristics

- High tensile strength.
- Excellent handling.
- Excellent knot security.
- Excellent tissue support.
- Smooth passage through the tissues.

Indications

Superficial soft tissue and/or ligation.

especially in

- General surgery
- Skin closure
- Gastrointestinal surgery
- Plastic surgery
- Gynaecology
- Urology


Contraindication

Cardiovascular surgery, ophthalmic surgery and neurosurgery.
Since it is an absorbable suture; It should not to use when an extended wound support is required.

Packaging

SURCAPRONE is available in USP sizes 6/0 through 2 (metric sizes 0.7-5.0).
SURCAPRONE is available with needles in one dozen boxes.
Three dozen suture boxes available on customer request.

TENSILE STRENGTH



Absorbable Polyglecaprone 25 Sutures

Round bodied 1/2 circle



Metric	USP	Length	Color	Order number
--------	-----	--------	-------	--------------

1	6/0	75 cm	violet	DR6016075
---	-----	-------	--------	-----------

1.5	5/0	75 cm	violet	DR5016075
-----	-----	-------	--------	-----------




1.5	4/0	75 cm	violet	DR4019075
-----	-----	-------	--------	-----------

2	3/0	75 cm	violet	DR3019075
---	-----	-------	--------	-----------



2	3/0	75 cm	violet	DR3025075
---	-----	-------	--------	-----------

3	2/0	75 cm	violet	DR2025075
---	-----	-------	--------	-----------



2	3/0	75 cm	violet	DR3030075
---	-----	-------	--------	-----------

3	2/0	75 cm	violet	DR2030075
---	-----	-------	--------	-----------


3.5	0	75 cm	violet	DR0030075
-----	---	-------	--------	-----------



3	2/0	75 cm	violet	DR2035075
---	-----	-------	--------	-----------

3.5	0	75 cm	violet	DR0035075
-----	---	-------	--------	-----------


4	1	75 cm	violet	DR0135075
---	---	-------	--------	-----------



3	2/0	75 cm	violet	DR2040075
---	-----	-------	--------	-----------

3.5	0	75 cm	violet	DR0040075
-----	---	-------	--------	-----------

4	1	90 cm	violet	DR0140090
---	---	-------	--------	-----------



3.5	0	75 cm	violet	DR0050075
-----	---	-------	--------	-----------

4	1	90 cm	violet	DR0150090
---	---	-------	--------	-----------

5	2	90 cm	violet	DR0250090
---	---	-------	--------	-----------

Reverse cutting 3/8 circle



Metric	USP	Length	Color	Order number
--------	-----	--------	-------	--------------

0.7	6/0	75 cm	violet	DC6016075
-----	-----	-------	--------	-----------

1	5/0	75 cm	violet	DC5016075
---	-----	-------	--------	-----------



1.5	4/0	75 cm	violet	DC4019075
-----	-----	-------	--------	-----------

2	3/0	75 cm	violet	DC3019075
---	-----	-------	--------	-----------



2	3/0	75 cm	violet	DC3025075
---	-----	-------	--------	-----------

3	2/0	75 cm	violet	DC2025075
---	-----	-------	--------	-----------



2	3/0	75 cm	violet	DC3030075
---	-----	-------	--------	-----------

3	2/0	75 cm	violet	DC2030075
---	-----	-------	--------	-----------

3.5	0	75 cm	violet	DC0030075
-----	---	-------	--------	-----------



3.5	0	75 cm	violet	DC0037075
-----	---	-------	--------	-----------


4	1	90 cm	violet	DC0137090
---	---	-------	--------	-----------

SURCAPRONE

Absorbable Polyglecaprone 25 Sutures

Reverse cutting 1/2 circle ▼	Metric	USP	Length	Color	Order number
	3	2/0	75 cm	violet	DC2035075
	3.5	0	75 cm	violet	DC0035075
	4	1	75 cm	violet	DC0135075
	3.5	0	90 cm	violet	DC0045090
	4	1	90 cm	violet	DC0145090
	5	2	90 cm	violet	DC0245090
Cutting straight ▲	Metric	USP	Length	Color	Order number
	2	3/0	75 cm	violet	DS3060075
	3	2/0	75 cm	violet	DS2060075
Heavy tapercut 5/8 circle ⊙	Metric	USP	Length	Color	Order number
	2	3/0	75 cm	violet	DT3025075
	2	3/0	75 cm	violet	DT3025075D
Without needle ○	Metric	USP	Length	Color	Order number
	2	3/0	150 cm	violet	DNN30150
	3	2/0	150 cm	violet	DNN20150
	3.5	0	150 cm	violet	DNN00150
	4	1	150 cm	violet	DNN01150
	5	2	150 cm	violet	DNN02150

Absorbable Polydioxanone Sutures



Type
Monofilament


Composition
Polydioxanone, a polymer made from polyester poly (p-dioxanone)

Coating
None

Color
Violet

Tissue reaction
Minimal

TENSILE STRENGTH



Absorption

The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. The material retains approximately 70% of its initial tensile strength at the end of fourth week. Approximately 50% of tensile strength remain after six weeks.

Distinctive characteristics

- High tensile strength.
- Excellent handling.
- Excellent knot security.
- Excellent tissue support.
- Smooth passage through the tissues.

Indications

- General soft tissue closing.
- especially in*
- Abdominal wall closure
- General surgery
- Gastrointestinal surgery
- Plastic surgery
- Ophthalmic surgery
- Gynaecology
- Orthopaedics

Contraindication















Cardiovascular surgery, neurosurgery and prosthetic devices.
Since it is an absorbable suture; It should not to use when an extended wound support is required.

Packaging

SURCANPONE is available in USP sizes 6/0 through 2 (metric sizes 0.7-5.0).
SURCANPONE is available with needles in one dozen boxes.
Three dozen suture boxes available on customer request.

SURCANONE

Absorbable Polydioxanone Sutures

Round bodied 1/2 circle		Metric	USP	Length	Color	Order number
	16MM	1	5/0	75 cm	violet	CR5016075
		1.5	4/0	75 cm	violet	CR4016075
	19MM	1.5	4/0	75 cm	violet	CR4019075
		2	3/0	75 cm	violet	CR3019075
	25MM	2	3/0	75 cm	violet	CR3025075
		3	2/0	75 cm	violet	CR2025075
	30MM	2	3/0	75 cm	violet	CR3030075
		3	2/0	75 cm	violet	CR2030075
		3.5	0	75 cm	violet	CR0030075
	35MM	3	2/0	75 cm	violet	CR2035075
		3.5	0	75 cm	violet	CR0035075
		4	1	75 cm	violet	CR0135075
	40MM	3	2/0	75 cm	violet	CR2040075
		3.5	0	75 cm	violet	CR0040075
		4	1	90 cm	violet	CR0140090
	50MM	3.5	0	75 cm	violet	CR0050075
		4	1	90 cm	violet	CR0150090
		5	2	90 cm	violet	CR0250090
Reverse cutting 3/8 circle		Metric	USP	Length	Color	Order number
	16MM	0.7	6/0	75 cm	violet	CC6016075
		1	5/0	75 cm	violet	CC5016075
	19MM	1.5	4/0	75 cm	violet	CC4019075
		2	3/0	75 cm	violet	CC3019075
	25MM	2	3/0	75 cm	violet	CC3025075
		3	2/0	75 cm	violet	CC2025075
	30MM	2	3/0	75 cm	violet	CC3030075
		3	2/0	75 cm	violet	CC2030075
		3.5	0	75 cm	violet	CC0030075
	40MM	3.5	0	75 cm	violet	CC0040075
		4	1	75 cm	violet	CC0140075


Absorbable Polydioxanone Sutures

Reverse cutting 1/2 circle ▼




Metric	USP	Length	Color	Order number
3	2/0	75 cm	violet	CC2037075
3.5	0	75 cm	violet	CC0037075
4	1	75 cm	violet	CC0137075


Round bodied heavy loop 1/2 circle ●



Metric	USP	Length	Color	Order number
3.5	0	150 cm	violet	CR0040150L
4	1	150 cm	violet	CR0140150L



Metric	USP	Length	Color	Order number
3.5	0	150 cm	violet	CR0045150L
4	1	150 cm	violet	CR0145150L



Metric	USP	Length	Color	Order number
3.5	0	150 cm	violet	CR0050150L
4	1	150 cm	violet	CR0150150L

Without needle ○



Metric	USP	Length	Color	Order number
4	1	150 cm	violet	CNN01150
5	2	150 cm	violet	CNN02150

SURSILK

Non Absorbable Silk Sutures

Type

Braided and coated multifilament

Composition

Braided fibres from the cocoon of the silkworm

Coating

Wax

Color

Black

Tissue reaction

Moderate

Absorption

SURSILK elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues.

Distinctive characteristics

High tensile strength.
High flexibility.
Excellent handling.
Excellent knot security.
Excellent tissue support for up to three months.
Smooth passage through the tissues.

Indications

General soft tissue closing and/or ligation.

especially in

General surgery
Skin closure
Gastrointestinal surgery
Oral surgery
Gynaecology
Ophthalmic surgery
Obstetrics

Contraindication

Since it is non absorbable suture; it is NOT recommended in urological procedures nor biliary tracts as this may lead to formation of stones around the suture.
Known allergies or sensitivities to silk.

Packaging

SURSILK is available in USP sizes 6/0 through 1 (metric sizes 0.7-4.0).
SURSILK is available with needles in one dozen boxes.
Three dozen suture boxes available on customer request.



Non Absorbable Silk Sutures

Round bodied 1/2 circle



Metric	USP	Length	Color	Order number
--------	-----	--------	-------	--------------

0.7	6/0	75 cm	black	SR6016075
-----	-----	-------	-------	-----------

1	5/0	75 cm	black	SR5016075
---	-----	-------	-------	-----------




1	5/0	75 cm	black	SR5019075
---	-----	-------	-------	-----------

1.5	4/0	75 cm	black	SR4019075
-----	-----	-------	-------	-----------



1.5	4/0	75 cm	black	SR4025075
-----	-----	-------	-------	-----------

2	3/0	75 cm	black	SR3025075
---	-----	-------	-------	-----------



2	3/0	75 cm	black	SR3030075
---	-----	-------	-------	-----------

3	2/0	75 cm	black	SR2030075
---	-----	-------	-------	-----------

3.5	0	75 cm	black	SR0030075
-----	---	-------	-------	-----------

4	1	75 cm	black	SR0130075
---	---	-------	-------	-----------



3	2/0	75 cm	black	SR2040075
---	-----	-------	-------	-----------

3.5	0	75 cm	black	SR0040075
-----	---	-------	-------	-----------

4	1	75 cm	black	SR0140075
---	---	-------	-------	-----------

Round bodied 3/8 circle



Metric	USP	Length	Color	Order number
--------	-----	--------	-------	--------------

2	3/0	75 cm	black	SR3026075
---	-----	-------	-------	-----------

3	2/0	75 cm	black	SR2026075
---	-----	-------	-------	-----------



3	2/0	75 cm	black	SR2032075
---	-----	-------	-------	-----------

3.5	0	75 cm	black	SR0032075
-----	---	-------	-------	-----------

Reverse cutting 3/8 circle



Metric	USP	Length	Color	Order number
--------	-----	--------	-------	--------------

0.7	6/0	75 cm	black	SC6016075
-----	-----	-------	-------	-----------

1	5/0	75 cm	black	SC5016075
---	-----	-------	-------	-----------



1.5	4/0	75 cm	black	SC4019075
-----	-----	-------	-------	-----------

2	3/0	75 cm	black	SC3019075
---	-----	-------	-------	-----------



2	3/0	75 cm	black	SC3025075
---	-----	-------	-------	-----------

3	2/0	75 cm	black	SC2025075
---	-----	-------	-------	-----------







2	3/0	75 cm	black	SC3030075
---	-----	-------	-------	-----------

3	2/0	75 cm	black	SC2030075
---	-----	-------	-------	-----------




3.5	0	75 cm	black	SC0030075
-----	---	-------	-------	-----------

Non Absorbable Silk Sutures



Reverse cutting 3/8 circle ▼

	Metric	USP	Length	Color	Order number
	3	2/0	75 cm	black	SC2037075
	3.5	0	75 cm	black	SC0037075
	4	1	75 cm	black	SC0137075
	3	2/0	75 cm	black	SC2040075
	3.5	0	75 cm	black	SC0040075
	4	1	75 cm	black	SC0140075
	5	2	75 cm	black	SC0240075
	3	2/0	75 cm	black	SC2045075
	3.5	0	75 cm	black	SC0045075
	4	1	75 cm	black	SC0145075
	5	2	75 cm	black	SC0245075
	3	2/0	75 cm	black	SC2050075
	3.5	0	75 cm	black	SC0050075
	4	1	75 cm	black	SC0150075
	5	2	75 cm	black	SC0250075

Reverse cutting 1/2 circle ▼

	Metric	USP	Length	Color	Order number
	1.5	4/0	75 cm	black	SC4018075
	2	3/0	75 cm	black	SC3018075
	2	3/0	75 cm	black	SC3026075
	3	2/0	75 cm	black	SC2026075
	3	2/0	75 cm	black	SC2038075
	3.5	0	75 cm	black	SC0038075
	4	1	75 cm	black	SC0138075

Tapercut 1/2 circle ⊕

	Metric	USP	Length	Color	Order number
	1.5	4/0	75 cm	black	ST4018075
	2	3/0	75 cm	black	ST3018075
	2	3/0	75 cm	black	ST3025075
	3	2/0	75 cm	black	ST2025075

Non Absorbable Silk Sutures

Cutting straight



Metric	USP	Length	Color	Order number
2	2/0	75 cm	black	SS2060075
3.5	0	75 cm	black	SS0060075
4	1	75 cm	black	SS0160075
5	2	75 cm	black	SS0260075

Without needle



Metric	USP	Length	Color	Order number
2	3/0	150 cm	black	SNN30150
3	2/0	150 cm	black	SNN20150
3.5	0	150 cm	black	SNN00150
4	1	150 cm	black	SNN01150
5	2	150 cm	black	SNN02150
2	3/0	6 X 45 cm	black	SNN30645
3	2/0	6 X 45 cm	black	SNN20645
3.5	0	6 X 45 cm	black	SNN00645
4	1	6 X 45 cm	black	SNN01645
5	2	6 X 45 cm	black	SNN02645
2	3/0	17 X 45 cm	black	SNN301745
3	2/0	17 X 45 cm	black	SNN201745
3.5	0	17 X 45 cm	black	SNN001745
4	1	17 X 45 cm	black	SNN011745
5	2	17 X 45 cm	black	SNN021745
2	3/0	13 X 60 cm	black	SNN301360
3	2/0	13 X 60 cm	black	SNN201360
3.5	0	13 X 60 cm	black	SNN001360
4	1	13 X 60 cm	black	SNN011360
5	2	13 X 60 cm	black	SNN021360

SURLYEST

Non Absorbable Polyester Sutures

Type

Braided and coated multifilament

Composition

Polyester - a polymer of polyethylene terephthalate

Coating

Silicone

Color

Green

Tissue reaction

Minimal

Absorption

Non absorbable, gradually encapsulated by connective tissue

Distinctive characteristics

Minimal tissue reaction.
Ideal for use in cardiovascular surgery
High flexibility.
Excellent knot security.
Excellent and permanent tissue support.
Good handling proprieties.
Smooth passage through the tissues.
Superior tensile strength.

Indications

General soft tissue closing and/or ligation.

especially in

Cardiovascular surgery
General surgery
Gastrointestinal surgery
Orthopaedics
Obstetrics
Ophthalmic surgery
Skin closure

Packaging









SURLONE is available in USP sizes 6/0 through 2 (metric sizes 0.7-5.0).







SURLONE is available with needles in one dozen boxes.

Three dozen suture boxes available on customer request.






Non Absorbable Polyester Sutures

Round bodied 1/2 circle		Metric	USP	Length	Color	Order number
	16MM	1.5	4/0	75 cm	green	LR4016075
		2	3/0	75 cm	green	LR3016075
	19MM	1.5	4/0	75 cm	green	LR4019075
		2	3/0	75 cm	green	LR3019075
	25MM	2	3/0	75 cm	green	LR3025075
		3	2/0	75 cm	green	LR2025075
	30MM	3	2/0	75 cm	green	LR2030075
		3.5	0	75 cm	green	LR0030075
		4	1	75 cm	green	LR0130075
	35MM	3	2/0	75 cm	green	LR2035075
		3.5	0	75 cm	green	LR0035075
		4	1	75 cm	green	LR0135075
	40MM	3.5	0	75 cm	green	LR0040075
		4	1	75 cm	green	LR0140075
		5	2	75 cm	green	LR0240075
	50MM	5	2	75 cm	green	LR0250075

Reverse cutting 3/8 circle		Metric	USP	Length	Color	Order number
	16MM	1	5/0	75 cm	green	LC5016075
		1.5	4/0	75 cm	green	LC4016075
	19MM	1.5	4/0	75 cm	green	LC4019075
		2	3/0	75 cm	green	LC3019075
	25MM	2	3/0	75 cm	green	LC3025075
		3	2/0	75 cm	green	LC2025075
	30MM	2	3/0	75 cm	green	LC3030075
		3	2/0	75 cm	green	LC2030075
		3.5	0	75 cm	green	LC0030075
	40MM	3.5	0	75 cm	green	LC0040075
		4	1	75 cm	green	LC0140075

SURLYEST

Non Absorbable Polyester Sutures

Reverse cutting 1/2 circle ▼	Metric	USP	Length	Color	Order number
	1.5	4/0	75 cm	green	LC4018075
	2	3/0	75 cm	green	LC3018075
	5	2	75 cm	green	LC0248075
Cutting straight ▲	Metric	USP	Length	Color	Order number
	2	3/0	75 cm	green	LS3060075
	3	2/0	75 cm	green	LS2060075
Without needle ○	Metric	USP	Length	Color	Order number
	2	3/0	150 cm	green	LNN30150
	3	2/0	150 cm	green	LNN20150
	3.5	0	150 cm	green	LNN00150
	4	1	150 cm	green	LNN01150
	5	2	150 cm	green	LNN02150
	2	3/0	6 X 45 cm	green	LNN30645
	3	2/0	6 X 45 cm	green	LNN20645
	3.5	0	6 X 45 cm	green	LNN00645
	4	1	6 X 45 cm	green	LNN01645
	5	2	6 X 45 cm	green	LNN02645
	2	3/0	3 X 50 cm	green	LNN30350
	3	2/0	3 X 50 cm	green	LNN20350
	3.5	0	3 X 50 cm	green	LNN00350
	4	1	3 X 50 cm	green	LNN01350
	5	2	3 X 50 cm	green	LNN02350

Non Absorbable Polypropylene Sutures

Type

Monofilament

Composition

Polypropylene, a polymer of propylene

Coating

None

Color

Blue

Tissue reaction

Moderate

Absorption

Non absorbable

Distinctive characteristics

- Low tissue reaction.
- Ideal biocompatibility.
- Excellent knot security.
- Excellent and permanent tissue support.
- Good handling proprieties.
- Smooth passage through the tissues.
- Absence of capillarity (no support of bacterial growth).

Indications

General soft tissue closing and/or ligation.

especially in

- General surgery
- Cardiovascular surgery
- Neurosurgery
- Ophthalmic surgery
- Plastic Syrgery
- Skin closure
- Gastrointestinal surgery
- Orthopaedics
- Gynaecology
- Obstetrics
- Microsurgery
- Hernias








Packaging






- SURLENE is available in USP sizes 6/0 through 1 (metric sizes 0.7-4.0).
- SURLENE is available with needles in one dozen boxes.
- Three dozen suture boxes available on customer request.










SURLENE


Non Absorbable Polypropylene Sutures



Round bodied 1/2 circle		Metric	USP	Length	Color	Order number
		0.7	6/0	75 cm	blue	PR6016075
		1	5/0	75 cm	blue	PR5016075
		1	5/0	75 cm	blue	PR5019075
		1.5	4/0	75 cm	blue	PR4019075
		1.5	4/0	75 cm	blue	PR4025075
		2	3/0	75 cm	blue	PR3025075
		2	3/0	75 cm	blue	PR3030075
		3	2/0	75 cm	blue	PR2030075
		3.5	0	75 cm	blue	PR0030075
		4	1	75 cm	blue	PR0130075
		3	2/0	75 cm	blue	PR2035075
		3.5	0	75 cm	blue	PR0035075
		4	1	75 cm	blue	PR0135075
		4	1	90 cm	blue	PR0140090
		5	2	90 cm	blue	PR0240090

Round bodied 1/2 circle		Metric	USP	Length	Color	Order number
Double armed 		0.7	6/0	90 cm	blue	PR6016090D
		1	5/0	90 cm	blue	PR5016090D
		1	5/0	90 cm	blue	PR5019090D
		1.5	4/0	90 cm	blue	PR4019090D
		2	3/0	90 cm	blue	PR3019090D
		1.5	4/0	90 cm	blue	PR4025090D
		2	3/0	90 cm	blue	PR3025090D
		3	2/0	90 cm	blue	PR2025090D
		3	2/0	90 cm	blue	PR2030075D

Non Absorbable Polypropylene Sutures

Reverse cutting 3/8 circle ▼	Metric	USP	Length	Color	Order number
	0.7	6/0	45 cm	blue	PC6016075
	1	5/0	45 cm	blue	PC5016075
	1.5	4/0	75 cm	blue	PC4019075
	2	3/0	75 cm	blue	PC3019075
	2	3/0	75 cm	blue	PC3025075
	3	2/0	75 cm	blue	PC2025075
	2	3/0	75 cm	blue	PC3030075
	3	2/0	75 cm	blue	PC2030075
	3.5	0	75 cm	blue	PC0030075
	2	3/0	75 cm	blue	PC3035075
	3	2/0	75 cm	blue	PC2035075
	3.5	0	75 cm	blue	PC0035075
	3	2/0	75 cm	blue	PC2040075
	3.5	0	90 cm	blue	PC0040090
	4	1	90 cm	blue	PC0140090
	4	1	90 cm	blue	PC0145090

Cutting straight ▲	Metric	USP	Length	Color	Order number
	2	3/0	75 cm	blue	PS3060075
	3	2/0	75 cm	blue	PS2060075

Tapercut 1/2 circle ⊕	Metric	USP	Length	Color	Order number
Double armed 	1.5	4/0	90 cm	blue	PT4016090D
		2	3/0	90 cm	blue
3		2/0	90 cm	blue	PT2025090D

SURLON

Non Absorbable Polyamide Sutures

Type

Monofilament

Composition

Extrusion of polyamide 6.0 or 6.6

Coating

None

Color

Blue
Black

Tissue reaction

Minimal

Absorption

Non absorbable, gradually encapsulated by connective tissue.

The thread mass diminishes, approximately 10% a year by rupture of chemical links (hydrolytic action).

Distinctive characteristics

Minimal tissue reaction (upon implementation).

Excellent visibility.

Excellent knot security.

Excellent and permanent tissue support.

Good handling proprieties.

Smooth passage through the tissues.

Absence of capillarity (no support of bacterial growth).

Indications

General soft tissue closing and/or ligation.

especially in

Blue

General surgery

Plastic Syrgery

Skin closure

Gastrointestinal surgery

Orthopaedics

Gynaecology

Obstetrics

Black

Plastic surgery

Ophthalmic surgery

Packaging

SURLON is available in USP sizes 6/0 through 2 (metric sizes 0.7-5.0).

SURLON is available with needles in one dozen boxes.








Three dozen suture boxes available on customer request.



Non Absorbable Polyamide Sutures







Round bodied 1/2 circle



	Metric	USP	Length	Color	Order number
	1	5/0	75 cm	black	NR5016075
	1.5	4/0	75 cm	black	NR4016075
	1	5/0	75 cm	black	NR5019075
	1.5	4/0	75 cm	black	NR4019075
	1.5	4/0	75 cm	black	NR4025075
	2	3/0	75 cm	black	NR3025075
	2	3/0	75 cm	black	NR3030075
	3	2/0	75 cm	black	NR2030075
	3.5	0	75 cm	black	NR0030075
	4	1	75 cm	black	NR0130075
	3	2/0	75 cm	black	NR2035075
	3.5	0	75 cm	black	NR0035075
	4	1	75 cm	black	NR0135075
	3	2/0	75 cm	black	NR2040075
	3.5	0	75 cm	black	NR0040075
	4	1	75 cm	black	NR0140075
	5	2	75 cm	black	NR0240075






Reverse cutting 3/8 circle




	Metric	USP	Length	Color	Order number
	0.7	6/0	75 cm	black	NC6016075
	1	5/0	75 cm	black	NC5016075
	1.5	4/0	75 cm	black	NC4019075
	2	3/0	75 cm	black	NC3019075
	2	3/0	75 cm	black	NC3025075
	3	2/0	75 cm	black	NC2025075
	2	3/0	75 cm	black	NC3030075
	3	2/0	75 cm	black	NC2030075
	3.5	0	75 cm	black	NC0030075
	3	2/0	75 cm	black	NC2040075
	3.5	0	75 cm	black	NC0040075
	4	1	90 cm	black	NC0140090
	5	2	90 cm	black	NC0240090


SURLON

Non Absorbable Polyamide Sutures

Reverse cutting 1/2 circle ▼	Metric	USP	Length	Color	Order number
	0.7	6/0	75 cm	black	NC6018075
	1	5/0	75 cm	black	NC5018075
	1.5	4/0	75 cm	black	NC4026075
	2	3/0	75 cm	black	NC3026075
	3	2/0	75 cm	black	NC2026075
	3	2/0	75 cm	black	NC2033075
	3.5	0	75 cm	black	NC0033075
	4	1	75 cm	black	NC0133075
	3	2/0	75 cm	black	NC2037075
	3.5	0	75 cm	black	NC0037075
	4	1	75 cm	black	NC0137075
	3.5	0	75 cm	black	NC0042075
	4	1	75 cm	black	NC0142075
	5	2	75 cm	black	NC0242075

Cutting straight ▲	Metric	USP	Length	Color	Order number
	2	3/0	75 cm	black	NS3060075
	2	2/0	75 cm	black	NS2060075

Non-absorbable Polypropylene Mesh



Composition

100% extruded monofilament polypropylene..

Color

Transparent

Tissue reaction

Moderate

Absorption

Non absorbable

Distinctive characteristics

Soft and pliable: For maximum patient comfort.
Light weight construction: Less foreign mass and flexible scar tissue.
Large pores: For optimal healing and scar formation.
Flexible: Recovers to original shape for any anatomic placement.
Full transparency: To see the tissue below.
Resist bacterial colonization: Construction allows fluid flow – through.
Long term stability: No shrinkage.
Pre Cut Sizes: Can be readily customized.

Indications

Inguinal hernia
Incisional hernia
Reconstruction of chest wall

Contraindication


SURGIMESH PRO should not be implanted in contaminated and infected areas as well as in children during the growth phase


Packaging







SURGIMESH PRO is available sterile, undyed in several sizes .
SURGIMESH PRO is individually packaged in two sealed pouches.


SURGIMESH PRO

Non-absorbable Polypropylene Mesh

Rectangle	Size	Indication	Order number
	2CM X 20CM	Prolapse	R0220
	5CM X 20CM	Prolapse	R0520
	5CM X 10CM	Open & laparoscopic hernia repair	R0510
	6CM X 11CM	Open & laparoscopic hernia repair	R0611
	8CM X 12CM	Open & laparoscopic hernia repair	R0812
	7.5CM X 15CM	Open & laparoscopic hernia repair	R7515
	8CM X 15CM	Open & laparoscopic hernia repair	R0815
	10CM X 13CM	Open & laparoscopic hernia repair	R1013
	10CM X 15CM	Open & laparoscopic hernia repair	R1015
	22.5CM X 35CM	Open & laparoscopic hernia repair	R22535
	25CM X 35CM	Open & laparoscopic hernia repair	R2535

Square	Size	Indication	Order number
	10CM X 10CM	Open & laparoscopic hernia repair	S1010
	15CM X 15CM	Open & laparoscopic hernia repair	S1515
	20CM X 20CM	Open & laparoscopic hernia repair	S2020
	30CM X 30CM	Open & laparoscopic hernia repair	S3030

Pre-Cut	Size	Indication	Order number
	5CM X 5CM	Indirect inguinal hernia	P0505
	7CM X 7CM	Indirect inguinal hernia	P0707
	5CM X 5CM	Indirect inguinal hernia	P0505H
	7CM X 7CM	Indirect inguinal hernia	P0707H
	4.5CM X 10CM	Open inguinal hernia repair	P4510
	6CM X 11CM	Open inguinal hernia repair	P0611
	4.5CM X 10CM	Open inguinal hernia repair	P4510H
	6CM X 11CM	Open inguinal hernia repair	P0611H
	7.5CM X 12.5CM	Rectum surgery	P75125
	8.5CM X 15CM	Rectum surgery	P8515
	12CM X 12CM	Hernia repair	P1212
	15CM X 15CM	Hernia repair	P1515

Plug	Size	Indication	Order number
	2.8CM X 3CM	Groin hernia surgical repair with anterior access	PL2803
	3.1CM X 3.4CM	Groin hernia surgical repair with anterior access	PL3134

Please feel free to contact us for other special sizes

Warranty of Quality ...

SURGIMEDIC offers high quality products to the market by making all processes from raw material to final product delivery without sacrificing quality.

It has adopted a continuous quality improvement policy that seeks good quality relations with customers before and after-sales, to receive feedbacks from customers and to improve quality in accordance with these feedbacks and legal requirements.

The manufacturing site has ISO class 8 cleanrooms with experienced staff that manufactures world quality surgical sutures and meshes.


With our ISO 13485:2016 and CE certificates, we guarantee the quality of our products.



IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.
© 2020 SURGIMEDIC. All rights reserved. SURGIMEDIC, SURGIMEDIC logo and Further, Together are trademarks of SURGIMEDIC.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products.

Reproduction in any form partial or otherwise is not permitted.



SURGIMEDIC

Quality. Commitment. Innovation.

Tel: **+20 (0) 12 22194006** | 104 - 106 A6 Industrial Zone, 10th of Ramadan 44916, Egypt
Web: <https://www.surgimedic.com> | Email: contact@surgimedic.com

